



Push[®] MetaGrip[®]

The Sophisticated Solution to a Common Problem

The most effective streamlined and durable brace for the relief of thumb CMC (basal joint) osteoarthritis.

Custom-fit the brace by squeezing the imbedded metal insert to snugly fit *your* thumb.

Wear the Push[®] MetaGrip[®] while cooking, golfing, gardening, playing tennis, driving, knitting, and during all other work/household/sports activities.

Advantages:

- Delivers long term durability, resisting abrasion
- Covers minimal surface of the palm
- Can be worn under a glove
- Can be worn underwater and washed in a washing machine
- Will not deform if left in a hot environment
- Is latex free

Indications:


- CMC Osteoarthritis



Measuring/Sizing

- To determine the appropriate Push MetaGrip size, measure the circumference of the palm just below the knuckles, excluding the thumb (see below). This circumferential measurement indicates hand size, but not thumb size. Please take this into consideration when choosing the appropriate size. If the measurement is near the size break point, only trying on both sizes can allow you to determine the best fit.
- **Ideal fitting around the thumb is accomplished by trying on the Push MetaGrip to assure proper fit. Professional measurement and fitting by a therapist/physician is suggested for optimal results.**

Size	Inches	Centimeters
1	6 3/8" - 7 3/4"	16 - 19.5
2	7 3/4" - 9"	19.5 - 22.5
3	9" - 10 1/4"	22.5 - 26



The Secret is in the Squeeze!

For the Push MetaGrip to be a successful solution to thumb CMC joint pain, it must be correctly SIZED and SQUEEZED to fit! Please visit our website www.HandLab.com for instructional videos on sizing, fitting, and caring for your Push MetaGrip.



HandLab

www.HandLab.com

support@handlab.com

Phone: 888-235-8221